S-SCORT[®] VX2

ENGLISH

PP 3-13

FRANÇAIS

PP 15-25

ITALIANO

PP 27-37

ESPAÑOL

PP 39-49

DEUTSCH

PP 51-61

**OPERATING INSTRUCTIONS
AND
MAINTENANCE MANUAL**

S-SCORT[®] VX2

Models

2310, 2310B, 2310V, 2310BV

2310-230, 2310B-230, 2310V-230, 2310BV-230

NEW GENERATION S-SCORT SUCTION FOR PREHOSPITAL CARE

30 LPM PERFORMANCE

DC INTERNAL POWER DC VEHICLE POWER

CONTROLLED BATTERY MAINTENANCE

REGULATED POWER TO SPARE

SSCOR, Inc.

11064 Randall St.

Sun Valley, CA 91352

USA

Tel +1 818 504 4054

Fax +1 818 504 6032

USA 800 434 5211

www.sscor.com

email: marketing@sscor.com

Table of Contents

S-SCORT VX2

Caution-Notice	4
General Description and Symbols	5
Getting Acquainted- Basic Instructions	6 – 7
General Specifications-Trouble Shooting	8
Warranty-Maintenance	9
Internal Access-Internal Components	10
Retention Bracket Mounting Instructions	9
Statement of Conformity	12

Caution-Notice

S-SCORT VX2

1. Federal law restricts this device to sale, distribution, and use by, or on the order of a physician, emergency medical technician, or other medical practitioner. For use by medical personnel trained in suctioning techniques and in the use of medical suction equipment. To be used only as directed by local EMS authorities.
2. This manual is restricted to the discussion of the use and maintenance of this device. It does not attempt to discuss professional techniques in suctioning procedures.
3. Operator should be thoroughly familiar with these operating instructions before this device is used.
4. Use a thumb control on all suction tips, especially when suctioning children and intubated patients.
5. The battery is protected from a deep discharge condition by shutting down after approximately 45 minutes ($\pm 10\%$) running time.
6. Do not use in the presence of flammable agents or anesthetics.
7. The shut-off valve in the canister lid will close down when the canister is full.
8. Install a new canister before testing for vacuum over 300mmHg to minimize the possibility of implosion, which can occur when a canister is aged or damaged.
9. All S-SCORT pumps must be reconnected to the charging source after each use and remain on charge until needed.

General Description

S-SCORT VX2


S-SCORT is a portable, 12V DC battery operated, suction device for resuscitation emergencies in Pre-Hospital Care.

CAUTION: S-SCORT suction units are designed for hospital crash cart suction, patient transport and emergency medical service. They are not designed or intended for use in extended procedures that require prolonged high vacuum/low airflow applications, as is the case in wound drainage or endoscopic use or in any other procedure that produces high vacuum levels within an occluded system for an extended period of time. Turn the S-SCORT suction unit off when it is not in use.

The unit receives power from its internal 12V DC battery or from the 12V system of the vehicle. The battery is recharged by an internal, dual mode battery charger which is activated by the AC-DC converter or direct connection to the vehicle electrical system. The charger monitors the battery, charges the battery only when necessary, shuts down the unit if battery is low and signals battery condition. A fully charged battery will power the unit for 45 minutes ($\pm 10\%$).

Description of Symbols

S-SCORT VX2

SYMBOL	LOCATION	MEANING
	Control panel	Pump is on
	Serialized label	Type B equipment
	Control panel	Battery level of charge Green LED = Full Green/Yellow LED = High Yellow led = Half Yellow/Red LED = Low Red LED = Battery depleted Put back on Charge
	Control panel	Push On / Push Off

S-SCORT® VX2 Models 2310, 2310B, 2310V, 2310BV, 2310-230, 2310B-230, 2310V-230, 2310BV-230



US Patent Pending

ETL listed, conforms to UL STD 2601-1, IEC 60601-1, EN 60601-1-2

©2003 SSCOR, Inc., S-SCORT, Duet, HI-D and Big Stick are registered trademarks of SSCOR, Inc.

Getting Acquainted

S-SCORT VX2

Getting Acquainted

S-SCORT VX2

Electrical operation: Constant battery charging is required.

DC BATTERY CHARGING FROM VEHICLE

Hard wire the charging cable to the hot DC circuit (in front of the master switch). Connect the electrical line cord to the S-SCORT by securely attaching the charging plug (1.1A) into the receptacle (1A). The vehicle electrical system furnishes power to the internal dual mode charger. The pump is powered by the vehicle current.

DC BATTERY CHARGING FROM AC POWER USING AC-DC CONVERTER

The 9Volt 1Amp AC Converter will not run the pump. If battery charging from AC power is required, connect the SSCOR AC-DC converter (1) to an AC outlet. The converter connects to the unit through the receptacle (1A) on the front panel. Do not attempt to run the pump with the AC converter attached to the pump. **Remove the AC converter connection from the pump prior to any suction procedure.**

BATTERY CHARGING VERIFICATION

Once the unit is connected to the charging source, check the control panel (2). The “power source connected” orange LED (3) shows a good connection and battery is charging.

TEST FOR VACUUM

- Press the ON/OFF switch (5) on the control panel to start the vacuum pump.
- Attach patient tubing to patient port (12) on collection canister (13).

Two Position Regulator Models

- Occlude patient tubing and push vacuum regulator (7) to the stop (full vacuum). Observe vacuum depressing canister lid. Pull vacuum regulator out to the stop to reduce vacuum to 125mmHg.

Variable Regulator Models

- Occlude patient tube and turn vacuum regulator (7) clockwise to the stop.
- Observe gauge (located above the vacuum regulator). The vacuum reading should rise to 300mmHg from the zero setting in 2 seconds. It should max out at approximately 525 mmHg.

NOTE 1: After the vacuum test, readjust the vacuum reading to your desired level.

NOTE 2: If the vacuum does not exceed 525mmHg, check for a leak in the system i.e. tubing, barb connections or loose canister lid (10,11,12,13). Refer internal vacuum problems to qualified personnel. After your test, reconnect the unit to the charging source; check the power light (3).

Running on DC Power

Turn the unit on and disconnect it from the charging source to test the function of the battery.

1. Check battery condition indicators (4), See page 5.

SPECIAL NOTICE: A blinking LED on the control panel indicates low battery capacity. This means the battery has been subjected to irreparable damage.

REPLACE BATTERY!

Keep the unit connected to the charging source at all times except for portable use.

General Specifications

S-SCORT VX2

CHARACTERISTICS	SPECIFICATIONS
Size	17"L x 9"H x 6"W
Weight	Approximately 9 lbs.
Vacuum Pump	12V DC oil-less diaphragm, 3.0A Exceeds 525 mmHg Clinical Airflow >= 30LPM
Regulator (7)	Controls negative pressure
Gauge (if applicable)	Calibrated in mmHg. Color-coded
Power Source: Battery (DC Power)	Rechargeable Sealed Lead Acid 45 (±10%) minute running time
Switch (5)	Off/Standby/On membrane switch
Collection Canister (13)	1200cc SSCOR part #48041
Patient tubing	Vinyl tubing 9/32"ID, 72"L SSCOR part #43200
Suction Tip	HI-D® . "Big Stick®". SSCOR part #44241

Trouble Shooting

S-SCORT VX2

PROBLEM	POSSIBLE CAUSE	CORRECTIVE ACTION
Does not function when switch is on (DC Power)	Battery discharged	Reconnect to charging source to activate charger
	Loose connections	Tighten connections
Blinking LED	Low battery capacity	Replace battery
Power source connected indicator light not lit	Power cord disconnected	Reconnect cord
	Converter defective	Check output of converter
No suction when pump is running	Vacuum line (10) loose	Tighten
	Canister (13) defective	Replace
	Canister lid (13) loose	Tighten
	Catheter thumb vent is open	Occlude
	Debris has collected in pump	Refer to maintenance section of this manual

Warranty

S-SCORT VX2

S-SCOR warrants that each new product is free from defects in material and workmanship under normal use and service for a period of one year from date of purchase.

If returned to SSCOR, we will arrange for repairs or replacement within the terms of the warranty. The defective product should be decontaminated and returned properly packaged and postage prepaid. Loss or damage in transit to the factory shall be at the purchaser's risk.

Please call 800/434-5211 for return authorization. Loss or damage in return shipment from SSCOR shall be at the purchaser's risk.

The warranty shall not apply to any SSCOR product which has been repaired by anyone other than an authorized SSCOR representative, or altered in any way so as, in SSCOR's judgment, to affect its safety or efficacy, nor which has been subject to misuse, negligence, or accident, nor which has had the serial number altered, effaced or removed.

Neither shall this warranty apply to any SSCOR product which has been connected otherwise than in accordance with the instructions furnished by SSCOR.

This warranty is in lieu of all other warranties expressed or implied and of all other obligations or liabilities on SSCOR's part, and SSCOR neither assumes, nor authorizes any representative or other persons to assume for it, any other liability in connection with the sale of SSCOR products.

Batteries and disposable collection canisters are excluded from this warranty.

S-SCORT VX2 with Charging Bracket Battery Replacement Policy

If the S-SCORT VX2 is purchased with the charging bracket SSCOR will replace the battery if the battery fails to operate the pump to specifications for a period of three years from the date of purchase. Verify that the pump runs to the specifications listed on page 7 of this Operations Manual.

Maintenance

S-SCORT VX2


- The unit is designed to suction contaminated fluids, which should be removed from the system immediately after the procedure. Dispose of the collection canister, patient tube and catheter and replace. If the procedure produced an excessive quantity of fluids, check the vacuum line (10) for evidence of moisture. Replace vacuum line if moisture is evident. (Use 3/8" ID tubing.)
- If the vacuum line between the pump and canister is moist, it is possible that fluids have reached the vacuum pump. Your engineering department will have to open the unit to check the condition of the pump.
- The capacity of the battery should be checked periodically by running the unit on battery only. Detach from the charging source and run the unit for 15 minutes and look for a blinking yellow light (page 7).
- If the Internal Dual Mode Charger appears defective, return the unit to the factory for repair. Do not attempt to repair the Internal Dual Mode Charger.
- Cleaning: Disconnect from power source prior to cleaning. Use a mild detergent. Do not use alcohol, solvents or cleaning solutions. For interior cleaning, use dry, compressed air.
- For technical assistance, call (800) 434-5211

Internal Access for Electrical and Pneumatics

S- SCORT VX2

WARNING: DO NOT ATTEMPT TO SERVICE THIS EQUIPMENT IF YOU ARE NOT A BIOMEDICAL TECHNICIAN

To Remove Canister Holder:

TO REMOVE CANISTER HOLDER:

Remove one 4" screw and two set screws on the canister bracket. Lift canister holder up and out of the main chassis.

OPEN/ CLOSE

Clamshell design enables access to components. All 6- 32 phillips head fastening screws located in back of the unit.

REMOVE MAIN CHASSIS SCREWS

- (A) 1 screw under the canister holder.
- (B) 2 screws in the back of the handle.
- (B) 3 screws in the back of the main chassis.
- (C) It should not be necessary to open the exhaust door unless fluids have entered the system. Slowly pull the unit apart at the handle to observe internal components.
- (E) PC BOARD- ELECTRICAL CIRCUITS
Do not attempt to service the PC Board. Return to SSCOR for service.
- (I) Battery Bracket.
- (J) Battery, 12V, Sealed Lead Acid.
- (K) Wiring connection to the control panel.
- (L) Vacuum pump. Return to SSCOR for service.
- (M) Exhaust barb
- (N) Vacuum barb

Charging/Retention Bracket Mounting Instructions

S-SCORT VX2

The Model 8310 Charging/Retention bracket can be mounted by the base or mounted vertically in order to hold the S-SCORT VX2 portable suction unit. Mount the unit so the indicator lights on the control panel are visible to the user. The unit is shipped ready to be mounted vertically, but if that does not fit your requirements, you can easily modify the unit to mount it by the base.

Note: the mounting bracket will hold the VX-2 in place to a force of ten 'G'. Select and prepare a mounting surface, and take care in the mounting procedure, so that the bracket and the pump will be secure at a force equal to, or greater than, 10 G

Vertical Mount

Remove the four 1/4 -20 x 1-1/2 hex head screws from hardware packet. Use the rectangle with the broken lines on the enclosed Mounting Template to mark the hole placement for a vertical mount. Use a 5/16" drill. Coat the screws with Loctite (not supplied) before securing the Charging/Retention Bracket. Put the four 1-1/2" screws through the holes in the charging bracket (Model 8310). Reattach the washers and nuts to the screws to secure the Charging/Retention bracket in place.

Base Mount (part #8314-7)

Use the rectangle with the solid lines on the enclosed Mounting Template to mark the hole placement for a base mount. Use a 5/16" drill. Remove the hardware from the hardware packet. Coat the screws with Loctite (not supplied) and place the screws through the bottom of the Bracket, Retention - Counter Mount (#8314-7) and through the holes you have drilled. Reattach the washers and nuts to secure the Bracket, Retention - Counter Mount (#8314-7) in place. Mount the charging bracket (Model 8310) to the base mount (#8314-7) using the hardware in the hardware packet. Coat the screws with Loctite (not supplied) before securing the Charging/Retention Bracket to the base mount bracket (#8314-7).

Electrical

Hard wire the DC Charging Cable to the DC Electrical System of the vehicle in front of the master switch. Plug the cable into the receptacle on the side of the bracket.

The ribbed wire of the charging cable is positive. The inside of the barrel connector to the charging bracket is positive. The front charging contact, closest to the canister holder, is positive and it matches the left charging contact on the charging bracket, which is positive.

DECLARATION OF CONFORMITY


*SSCOR, Incorporated
11064 Randall Street
Sun Valley, CA 91352
United States of America
Tel: (818) 504-4054
Fax: (818) 504-6032*

Type of Equipment: Portable Aspiration Pump
Model Name: S-SCORT VX2
Model Type Number: 2310, 2310B, 2310V, 2310BV, 2310-230,
2310B-230, 2310V-230, 2310BV-230

We herewith declare that the above mentioned product meets with the following requirements of the EC Directive concerning Medical Devices – 93/42/EEC

Safety: EN 60601-1: 1990
With A1 & A12:1993, A2:1995 & A13:1996
EMC: EN 60601-1-2: 1993
Classification Class IIa
Rules of Classification Annex IX, Rules 6 & 11
Annex used for CE Declaration Annex VI & Annex VII
Notified Body SEMKO, AB

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Gemeldete Behörde SEMKO, AB
**Dynamed UK Ltd